

IFN Gamma / Interferon Gamma Rabbit anti-Guinea pig Polyclonal (FITC) Antibody - LS-C210760 - LSBio	
CatalogID:	LS-C210760
Target:	interferon, gamma (IFNG)
Synonyms:	IFNG Antibody, IFG Antibody, Immune interferon Antibody, IFI Antibody, Interferon gamma Antibody, IFN-gamma Antibody, Interferon, gamma Antibody
Family / Subfamily:	Interferon / not assigned-Interferon
Host	IFNG antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	IFN Gamma / Interferon Gamma antibody was raised against Guinea pig
Antigen Type:	Recombinant protein
Immunogen:	IFN Gamma / Interferon Gamma antibody was raised against recombinant Guinea Pig Interferon gamma protein.
Specificity:	Guinea pig IFN Gamma / Interferon Gamma
Reactivity:	Guinea pig
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Concentration:	0.25 mg/ml
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	