

TTR / Transthyretin Guinea pig anti-Human Polyclonal (FITC) Antibody - LS-C210822 - LSBio

CatalogID:	LS-C210822
Target:	transthyretin (TTR)
Synonyms:	TTR Antibody, ATTR Antibody, CTS1 Antibody, CTS Antibody, HsT2651 Antibody, PALB Antibody, Prealbumin Antibody, TBPA Antibody, Transthyretin Antibody, Carpal tunnel syndrome 1 Antibody, Prealbumin, amyloidosis type I Antibody, Thyroxine-binding prealbumin Antibody
Host	TTR antibody was produced in Guinea pig
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	TTR / Transthyretin antibody was raised against Human
Antigen Type:	Purified protein
Immunogen:	TTR / Transthyretin antibody was raised against human Prealbumin protein.
Specificity:	Human TTR / Transthyretin
Reactivity:	Human
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences