

IL1B Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C210838 - LSBio

CatalogID:	LS-C210838
Target:	interleukin 1, beta (IL1B)
Synonyms:	IL1B Antibody, IL-1B Antibody, IL-1 beta Antibody, IL1-BETA Antibody, Preinterleukin 1 beta Antibody, Pro-interleukin-1-beta Antibody, Catabolin Antibody, IL-1 Antibody, IL1F2 Antibody, Interleukin 1, beta Antibody, Interleukin-1 beta Antibody
Family / Subfamily:	IL-1 / not assigned-IL-1
Host	IL1B antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	IL1B antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	IL1B antibody was raised against recombinant human Interleukin-1 beta protein.
Specificity:	Human IL1B / IL-1 Beta
Reactivity:	Human
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences