

EIF4A2 Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C210939 - LSBio	
CatalogID:	LS-C210939
Target:	eukaryotic translation initiation factor 4A2 (EIF4A2)
Synonyms:	EIF4A2 Antibody, BM-010 Antibody, EIF-4A-II Antibody, EIF4A Antibody, EIF4A-II Antibody, EIF4F Antibody, DDX2B Antibody
Host	EIF4A2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	EIF4A2 antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	EIF4A2 antibody was raised against recombinant human Eukaryotic initiation factor 4A-II protein.
Specificity:	Human EIF4A2
Reactivity:	Human
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	