

PFDN5 / MM1 Rabbit anti-Human Polyclonal (Biotin) Antibody - LS-C211119 - LSBio

| | |
|-----------------------------|---|
| CatalogID: | LS-C211119 |
| Target: | prefoldin subunit 5 (PFDN5) |
| Synonyms: | PFDN5 Antibody, C-Myc-binding protein Mm-1 Antibody, C-myc binding protein Antibody, MM-1 Antibody, Myc modulator 1 Antibody, Myc modulator-1 Antibody, Prefoldin 5 Antibody, MM1 Antibody, PFD5 Antibody, Prefoldin subunit 5 Antibody |
| Host | PFDN5 antibody was produced in Rabbit |
| Clonality: | Polyclonal |
| Isotype: | IgG |
| Conjugations: | Biotin |
| Immunogen Species: | PFDN5 / MM1 antibody was raised against Human |
| Antigen Type: | Recombinant protein |
| Immunogen: | PFDN5 / MM1 antibody was raised against recombinant human Prefoldin subunit 5 protein. |
| Specificity: | Human PFDN5 / MM1 |
| Reactivity: | Human |
| Purification: | Caprylic acid and ammonium sulfate precipitation |
| Presentation: | PBS. |
| Recommended Storage: | Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles. |
| Uses: | IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher) |
| Size: | 100 µg |
| Requested From: | Japan |

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences