

CSNK1E / CK1 Epsilon Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C211142 - LSBio	
CatalogID:	LS-C211142
Target:	casein kinase 1, epsilon (CSNK1E)
Synonyms:	CSNK1E Antibody, Casein kinase 1 epsilon Antibody, Casein kinase i epsilon Antibody, CK1 epsilon Antibody, CKIE Antibody, CKIepsilon Antibody, CKI-epsilon Antibody, HCKIE Antibody, Casein kinase 1, epsilon Antibody
Family / Subfamily:	Protein Kinase / Casein kinase I
Host	CSNK1E antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CSNK1E / CK1 Epsilon antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	CSNK1E / CK1 Epsilon antibody was raised against recombinant human Casein kinase I isoform epsilon protein.
Specificity:	Human CSNK1E / CK1 Epsilon
Reactivity:	Human
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	