

LCN2 / Lipocalin 2 / NGAL Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C211186 - LSBio	
CatalogID:	LS-C211186
Target:	lipocalin 2 (LCN2)
Synonyms:	LCN2 Antibody, 24p3 Antibody, HNL Antibody, Lipocalin 2 Antibody, Lipocalin-2 Antibody, MSF1 Antibody, Siderocalin LCN2 Antibody, NGAL Antibody, p25 Antibody, Lipocalin 2 (oncogene 24p3) Antibody, Oncogene 24p3 Antibody
Host	LCN2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	LCN2 / Lipocalin 2 / NGAL antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	LCN2 / Lipocalin 2 / NGAL antibody was raised against recombinant human NGAL protein.
Specificity:	Human LCN2 / Lipocalin 2 / NGAL
Reactivity:	Human
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	