

NPPB / BNP Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C211189 - LSBio

CatalogID:	LS-C211189
Target:	natriuretic peptide B (NPPB)
Synonyms:	NPPB Antibody, BNP Antibody, Natriuretic peptides B Antibody, Natriuretic peptide B Antibody, Natriuretic protein Antibody, Brain type natriuretic peptide Antibody
Host	NPPB antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	NPPB / BNP antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	NPPB / BNP antibody was raised against recombinant human NT-proBNP protein.
Specificity:	Human NPPB / BNP
Reactivity:	Human
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences