

ANXA5 / Annexin V Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C211251 - LSBio

CatalogID:	LS-C211251
Target:	annexin A5 (ANXA5)
Synonyms:	ANXA5 Antibody, ANX5 Antibody, Anchorin CII Antibody, Calphobindin I Antibody, CBP-I Antibody, ENX2 Antibody, PP4 Antibody, RPRGL3 Antibody, Lipocortin V Antibody, Thromboplastin inhibitor Antibody, VAC-alpha Antibody, Vascular anticoagulant-alpha Antibody, PAP-I Antibody, Annexin A5 Antibody, Annexin V Antibody, Annexin-5 Antibody, Endonexin II Antibody
Host	ANXA5 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ANXA5 / Annexin V antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	ANXA5 / Annexin V antibody was raised against recombinant human Annexin A5 protein.
Specificity:	Human ANXA5 / Annexin V
Reactivity:	Human
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences