

TUBB3 / Tubulin Beta 3 Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C211293 - LSBio

CatalogID:	LS-C211293
Target:	tubulin, beta 3 class III (TUBB3)
Synonyms:	TUBB3 Antibody, CDCBM Antibody, Beta III Tubulin Antibody, Class III beta-tubulin Antibody, TUBB4 Antibody, Tubulin, beta 3 Antibody, Tubulin beta-III Antibody, Tubulin beta-3 chain Antibody, Tubulin beta-4 chain Antibody, Tubulin, beta 3 class III Antibody, CFEOM3A Antibody
Family / Subfamily:	Tubulin
Host	TUBB3 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	TUBB3 / Tubulin Beta 3 antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	TUBB3 / Tubulin Beta 3 antibody was raised against recombinant human Tubulin beta-3 chain protein.
Specificity:	Human TUBB3 / Tubulin Beta 3
Reactivity:	Human
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences