

CANT1 Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C211305 - LSBio	
CatalogID:	LS-C211305
Target:	calcium activated nucleotidase 1 (CANT1)
Synonyms:	CANT1 Antibody, Apyrase homolog Antibody, DBQD Antibody, SCAN-1 Antibody, SHAPY Antibody, SCAN1 Antibody
Host	CANT1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CANT1 antibody was raised against Human
Reactivity:	Human
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	