

ERAS Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C211340 - LSBio

| | |
|-----------------------------|--|
| CatalogID: | LS-C211340 |
| Target: | ES cell expressed Ras (ERAS) |
| Synonyms: | ERAS Antibody, ES cell expressed Ras Antibody, E-Ras Antibody, GTPase ERas Antibody, HRASP Antibody, HRAS2 Antibody, Small GTPase protein E-Ras Antibody |
| Host | ERAS antibody was produced in Rabbit |
| Clonality: | Polyclonal |
| Isotype: | IgG |
| Conjugations: | Fluorescein (FITC) |
| Immunogen Species: | ERAS antibody was raised against Human |
| Antigen Type: | Recombinant protein |
| Immunogen: | ERAS antibody was raised against recombinant human GTPase ERas protein. |
| Specificity: | Human ERAS |
| Reactivity: | Human |
| Purification: | Caprylic acid and ammonium sulfate precipitation |
| Presentation: | PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol. |
| Recommended Storage: | Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles. |
| Uses: | IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher) |
| Size: | 100 µg |
| Requested From: | Japan |

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences