

Bikunin / AMBP Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C211380 - LSBio

CatalogID:	LS-C211380
Target:	alpha-1-microglobulin/bikunin precursor (AMBP)
Synonyms:	AMBP Antibody, A1M Antibody, Alpha-1-microglobulin Antibody, Bikunin Antibody, HCP Antibody, IATIL Antibody, ITIL Antibody, HI30 Antibody, ITI Antibody, ITILC Antibody, Protein AMBP Antibody, Uronic-acid-rich protein Antibody, UTI Antibody, Trypstatin Antibody, Uristatin Antibody, EDC1 Antibody, Growth-inhibiting protein 19 Antibody, Protein HC Antibody
Host	AMBP antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Bikunin / AMBP antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	Bikunin / AMBP antibody was raised against recombinant human Protein AMBP protein.
Specificity:	Human Bikunin / AMBP
Reactivity:	Human
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences