

SERPINE1 / PAI-1 Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C211384 - LSBio

CatalogID:	LS-C211384
Target:	serpin peptidase inhibitor, clade E (nexin, plasminogen activator inhibitor type 1), member 1 (SERPINE1)
Synonyms:	SERPINE1 Antibody, PAI1 Antibody, PLANH1 Antibody, PAI-1 Antibody, Serpin E1 Antibody
Family / Subfamily:	Serine Protease Inhibitor
Host	SERPINE1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	SERPINE1 / PAI-1 antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	SERPINE1 / PAI-1 antibody was raised against recombinant human Plasminogen activator inhibitor 1 protein.
Specificity:	Human SERPINE1 / PAI-1
Reactivity:	Human
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences