

IL1RN Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C211411 - LSBio

CatalogID:	LS-C211411
Target:	interleukin 1 receptor antagonist (IL1RN)
Synonyms:	IL1RN Antibody, Anakinra Antibody, DIRA Antibody, IL-1ra Antibody, IL-1RN Antibody, IL1F3 Antibody, IL1RA Antibody, MVCD4 Antibody, IL1 inhibitor Antibody, IL1RN (IL1F3) Antibody, IRAP Antibody, ICIL-1RA Antibody, IL-1ra3 Antibody
Family / Subfamily:	IL-1 / not assigned-IL-1
Host	IL1RN antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	IL1RN antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	IL1RN antibody was raised against recombinant human Interleukin-1 receptor antagonist protein.
Specificity:	Human IL1RN / IL-1RA
Reactivity:	Human
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences