

EGF Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C211516 - LSBio

CatalogID:	LS-C211516
Target:	epidermal growth factor (EGF)
Synonyms:	EGF Antibody, Beta-urogastrone Antibody, HOMG4 Antibody, URG Antibody, Epidermal growth factor Antibody, Pro-epidermal growth factor Antibody
Host	EGF antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	EGF antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	EGF antibody was raised against recombinant human Pro-epidermal growth factor protein.
Specificity:	Human EGF
Reactivity:	Human
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences