

**TNNT2 / CTNT Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C211524 - LSBio**

<b>CatalogID:</b>	LS-C211524
<b>Target:</b>	troponin T type 2 (cardiac) (TNNT2)
<b>Synonyms:</b>	TNNT2 Antibody, Cardiomyopathy, hypertrophic 2 Antibody, CMPD2 Antibody, CMD1D Antibody, CMH2 Antibody, Cardiac muscle troponin T Antibody, Cardiac Troponin T Antibody, LVNC6 Antibody, TnTC Antibody, Troponin T type 2 (cardiac) Antibody, Troponin T, cardiac muscle Antibody, Troponin T2, cardiac Antibody, CTnT Antibody, RCM3 Antibody
<b>Host</b>	TNNT2 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	TNNT2 / CTNT antibody was raised against Human
<b>Antigen Type:</b>	Recombinant protein
<b>Immunogen:</b>	TNNT2 / CTNT antibody was raised against recombinant human Troponin T, cardiac muscle protein.
<b>Specificity:</b>	Human TNNT2 / CTnT
<b>Reactivity:</b>	Human
<b>Purification:</b>	Caprylic acid and ammonium sulfate precipitation
<b>Presentation:</b>	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
<b>Recommended Storage:</b>	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
<b>Uses:</b>	Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences