

THBS1 / Thrombospondin-1 Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C211537 - LSBio	
CatalogID:	LS-C211537
Target:	thrombospondin 1 (THBS1)
Synonyms:	THBS1 Antibody, Thrombospondin-1 Antibody, Thrombospondin-1p180 Antibody, TSP Antibody, TSP1 Antibody, THBS-1 Antibody, THBS Antibody, Thrombospondin 1 Antibody, TSP-1 Antibody
Host	THBS1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	THBS1 / Thrombospondin-1 antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	THBS1 / Thrombospondin-1 antibody was raised against recombinant human Thrombospondin-1 protein.
Specificity:	Human THBS1 / Thrombospondin-1
Reactivity:	Human
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	