


PHGDH / PGDH Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C211553 - LSBio

CatalogID:	LS-C211553
Target:	phosphoglycerate dehydrogenase (PHGDH)
Synonyms:	PHGDH Antibody, 3-PGDH Antibody, 3PGDH Antibody, PGDH Antibody, Phosphoglycerate dehydrogenase Antibody, PDG Antibody, PGAD Antibody, PGDH3 Antibody
Host	PHGDH antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	PHGDH / PGDH antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	PHGDH / PGDH antibody was raised against recombinant human D-3-phosphoglycerate dehydrogenase protein.
Specificity:	Human PHGDH / PGDH
Reactivity:	Human
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg

Western Blot Image:



Western blot of D-3-phosphoglycerate dehydrogenase antibody at 2 ug/ml. Lane 1: EC109 whole cell lysate. Lane 2: 293T whole cell lysate. Secondary: Goat polyclonal to Rabbit IgG at 1:15000 dilution. Predicted band size: 59 kDa. Observed band size: 70 kDa Additional bands at: 100 kDa. 36 kDa. We are unsure as to the identity of these extra bands.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences