

NELFB / COBRA1 Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C211612 - LSBio

CatalogID:	LS-C211612
Target:	negative elongation factor complex member B (NELFB)
Synonyms:	NELFB Antibody, Cofactor of BRCA1 Antibody, COBRA1 Antibody, RP13-122B23.3 Antibody, KIAA1182 Antibody, Negative elongation factor B Antibody, NELF-B Antibody
Host	NELFB antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	NELFB / COBRA1 antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	NELFB / COBRA1 antibody was raised against recombinant human Negative elongation factor B protein.
Specificity:	Human COBRA1
Reactivity:	Human
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences