

AP2S1 Rabbit anti-Human Polyclonal Antibody - LS-C211718 - LSBio

CatalogID:	LS-C211718
Target:	adaptor-related protein complex 2, sigma 1 subunit (AP2S1)
Synonyms:	AP2S1 Antibody, AP17-DELTA Antibody, AP-2 complex subunit sigma Antibody, AP17 Antibody, CLAPS2 Antibody, HA2 17 kDa subunit Antibody, Sigma2-adaptin Antibody
Host	AP2S1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Immunogen Species:	AP2S1 antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	AP2S1 antibody was raised against recombinant human AP-2 complex subunit sigma protein.
Specificity:	Human AP2S1
Reactivity:	Human
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences