

AP2S1 Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C211719 - LSBio	
<b>CatalogID:</b>	LS-C211719
<b>Target:</b>	adaptor-related protein complex 2, sigma 1 subunit (AP2S1)
<b>Synonyms:</b>	AP2S1 Antibody, AP17-DELTA Antibody, AP-2 complex subunit sigma Antibody, AP17 Antibody, CLAPS2 Antibody, HA2 17 kDa subunit Antibody, Sigma2-adaptin Antibody
<b>Host</b>	AP2S1 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	AP2S1 antibody was raised against Human
<b>Antigen Type:</b>	Recombinant protein
<b>Immunogen:</b>	AP2S1 antibody was raised against recombinant human AP-2 complex subunit sigma protein.
<b>Specificity:</b>	Human AP2S1
<b>Reactivity:</b>	Human
<b>Purification:</b>	Caprylic acid and ammonium sulfate precipitation
<b>Presentation:</b>	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
<b>Recommended Storage:</b>	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
<b>Uses:</b>	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b>	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	