

IL4 Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C211747 - LSBio	
CatalogID:	LS-C211747
Target:	interleukin 4 (IL4)
Synonyms:	IL4 Antibody, B cell growth factor 1 Antibody, B-cell stimulatory factor 1 Antibody, B_cell stimulatory factor 1 Antibody, BCGF1 Antibody, BSF-1 Antibody, Binetrakin Antibody, BCGF-1 Antibody, Interleukin 4 Antibody, Interleukin-4 Antibody, Pitrakinra Antibody, BSF1 Antibody, IL-4 Antibody
Family / Subfamily:	Hormone / not assigned-Hormone
Host	IL4 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	IL4 antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	IL4 antibody was raised against recombinant human Interleukin-4 protein.
Specificity:	Human IL4 / IL-4
Reactivity:	Human
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/27/2014 © 2014 LifeSpan BioSciences	