

LYPLA2 Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C211832 - LSBio

CatalogID:	LS-C211832
Target:	lysophospholipase II (LYPLA2)
Synonyms:	LYPLA2 Antibody, Acyl-protein thioesterase 2 Antibody, APT-2 Antibody, APT2 Antibody, DJ886K2.4 Antibody, LysoPLA II Antibody, LPL-II Antibody, Lysophospholipase II Antibody
Family / Subfamily:	Phospholipase / not assigned-Phospholipase
Host	LYPLA2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	LYPLA2 antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	LYPLA2 antibody was raised against recombinant human Acyl-protein thioesterase 2 protein.
Specificity:	Human LYPLA2
Reactivity:	Human
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences