

GSTP1 / GST Pi Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C211861 - LSBio	
CatalogID:	LS-C211861
Target:	glutathione S-transferase pi 1 (GSTP1)
Synonyms:	GSTP1 Antibody, Deafness, X-linked 7 Antibody, DFN7 Antibody, GSTP1-1 Antibody, Glutathione S-transferase P Antibody, GST Pi Antibody, GSTP Antibody, FAES3 Antibody, Glutathione S-transferase pi 1 Antibody, GST class-pi Antibody, GST3 Antibody, PI Antibody
Host	GSTP1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	GSTP1 / GST Pi antibody was raised against Human
Reactivity:	Human
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014 © 2014 LifeSpan BioSciences	