

IP-10 / CXCL10 Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C211950 - LSBio

CatalogID:	LS-C211950
Target:	chemokine (C-X-C motif) ligand 10 (CXCL10)
Synonyms:	CXCL10 Antibody, C-X-C motif chemokine 10 Antibody, Gamma-IP10 Antibody, IP-10 Antibody, Gamma IP10 Antibody, GIP-10 Antibody, Mob-1 Antibody, Small-inducible cytokine B10 Antibody, Crg-2 Antibody, IFI10 Antibody, INP10 Antibody, SCYB10 Antibody
Family / Subfamily:	Intercrine / Intercrine alpha
Host	CXCL10 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	IP-10 / CXCL10 antibody was raised against Human
Reactivity:	Human
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences