

Ferritin Heavy Chain / FTH1 Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C211953 - LSBio	
CatalogID:	LS-C211953
Target:	ferritin, heavy polypeptide 1 (FTH1)
Synonyms:	FTH1 Antibody, Apoferritin Antibody, Ferritin H subunit Antibody, Ferritin heavy chain Antibody, Ferritin, heavy polypeptide 1 Antibody, FHC Antibody, FTHL6 Antibody, PIG15 Antibody, PLIF Antibody, FTH Antibody
Host	FTH1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Ferritin Heavy Chain / FTH1 antibody was raised against Human
Reactivity:	Human
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/27/2014 © 2014 LifeSpan BioSciences	