

COL4A1 / Collagen IV Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C211979 - LSBio

CatalogID:	LS-C211979
Target:	collagen, type IV, alpha 1 (COL4A1)
Synonyms:	COL4A1 Antibody, Arresten Antibody, COL4A1 NC1 domain Antibody, Collagen alpha-1(IV) chain Antibody, Collagen, type IV, alpha 1 Antibody, HANAC Antibody, ICH Antibody, POREN1 Antibody
Family / Subfamily:	Collagen / not assigned-Collagen
Host	COL4A1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	COL4A1 / Collagen IV antibody was raised against Human
Reactivity:	Human
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences