

ORM1 / Orosomuroid Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C211997 - LSBio

CatalogID:	LS-C211997
Target:	orosomuroid 1 (ORM1)
Synonyms:	ORM1 Antibody, AGP-A Antibody, Alpha-1-acid glycoprotein 1 Antibody, AGP 1 Antibody, AGP1 Antibody, Alpha-1-acid glycoprotein 2 Antibody, OMD 1 Antibody, Orosomuroid Antibody, Orosomuroid-1 Antibody, ORM Antibody, Orosomuroid 1 Antibody
Host	ORM1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ORM1 / Orosomuroid antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	ORM1 / Orosomuroid antibody was raised against recombinant human Alpha-1-acid glycoprotein 1 protein.
Specificity:	Human ORM1 / Orosomuroid
Reactivity:	Human
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences