

IL6 Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C212007 - LSBio	
CatalogID:	LS-C212007
Target:	interleukin 6 (interferon, beta 2) (IL6)
Synonyms:	IL6 Antibody, B-cell stimulatory factor 2 Antibody, BSF-2 Antibody, BSF2 Antibody, CTL differentiation factor Antibody, Hybridoma growth factor Antibody, IFN-beta-2 Antibody, Interleukin-6 Antibody, Interferon beta-2 Antibody, Interleukin BSF-2 Antibody, IFNB2 Antibody, B-cell differentiation factor Antibody, CDF Antibody, HSF Antibody, IL-6 Antibody
Family / Subfamily:	Hormone / not assigned-Hormone
Host	IL6 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	IL6 antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	IL6 antibody was raised against recombinant human Interleukin-6 protein.
Specificity:	Human IL6 / IL-6
Reactivity:	Human
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	