

FN1 / Fibronectin Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C212044 - LSBio

CatalogID:	LS-C212044
Target:	fibronectin 1 (FN1)
Synonyms:	FN1 Antibody, CIG Antibody, Cold-insoluble globulin Antibody, Fibronectin 1 Antibody, FINC Antibody, LETS Antibody, GFND Antibody, ED-B Antibody, Fibronectin Antibody, FN Antibody, FNZ Antibody, GFND2 Antibody, Migration-stimulating factor Antibody, MSF Antibody
Host	FN1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	FN1 / Fibronectin antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	FN1 / Fibronectin antibody was raised against recombinant human Fibronectin protein.
Specificity:	Human FN1 / Fibronectin
Reactivity:	Human
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences