

FOXP3 Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C212082 - LSBio	
<b>CatalogID:</b>	LS-C212082
<b>Target:</b>	forkhead box P3 (FOXP3)
<b>Synonyms:</b>	FOXP3 Antibody, AIID Antibody, DIETER Antibody, Forkhead box protein P3 Antibody, JM2 Antibody, Forkhead box P3 Antibody, FOXP3delta7 Antibody, XPID Antibody, PIDX Antibody, SCURFIN Antibody
<b>Host</b>	FOXP3 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	FOXP3 antibody was raised against Human
<b>Antigen Type:</b>	Recombinant protein
<b>Immunogen:</b>	FOXP3 antibody was raised against recombinant human Forkhead box protein P3 protein.
<b>Specificity:</b>	Human FOXP3
<b>Reactivity:</b>	Human
<b>Purification:</b>	Caprylic acid and ammonium sulfate precipitation
<b>Presentation:</b>	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
<b>Recommended Storage:</b>	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
<b>Uses:</b>	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	