

KLK10 / Kallikrein 10 Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C212087 - LSBio	
CatalogID:	LS-C212087
Target:	kallikrein-related peptidase 10 (KLK10)
Synonyms:	KLK10 Antibody, NES1 Antibody, Kallikrein-10 Antibody, Protease, serine-like, 1 Antibody, PRSSL1 Antibody, Kallikrein 10 Antibody, Protease serine-like 1 Antibody
Family / Subfamily:	Protease / Serine S1
Host	KLK10 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	KLK10 / Kallikrein 10 antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	KLK10 / Kallikrein 10 antibody was raised against recombinant human Kallikrein-10 protein.
Specificity:	Human KLK10 / Kallikrein 10
Reactivity:	Human
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	