

**GM2A Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C212104 - LSBio**

<b>CatalogID:</b>	LS-C212104
<b>Target:</b>	GM2 ganglioside activator (GM2A)
<b>Synonyms:</b>	GM2A Antibody, Ganglioside GM2 activator Antibody, GM2 ganglioside activator Antibody, GM2-AP Antibody, SAP-3 Antibody
<b>Host</b>	GM2A antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	GM2A antibody was raised against Human
<b>Antigen Type:</b>	Recombinant protein
<b>Immunogen:</b>	GM2A antibody was raised against recombinant human Ganglioside GM2 activator protein.
<b>Specificity:</b>	Human GM2A
<b>Reactivity:</b>	Human, Mouse
<b>Purification:</b>	Caprylic acid and ammonium sulfate precipitation
<b>Presentation:</b>	PBS.
<b>Recommended Storage:</b>	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
<b>Uses:</b>	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences