

Adiponectin Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C212133 - LSBio

CatalogID:	LS-C212133
Target:	adiponectin, C1Q and collagen domain containing (ADIPOQ)
Synonyms:	ADIPOQ Antibody, Adiponectin Antibody, ADPN Antibody, ACDC Antibody, ADIPQTL1 Antibody, APM-1 Antibody, APM1 Antibody, ACRP30 Antibody, Gelatin-binding protein 28 Antibody, GBP28 Antibody, Gelatin-binding protein Antibody
Host	ADIPOQ antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Adiponectin antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	Adiponectin antibody was raised against recombinant human Adiponectin protein.
Specificity:	Human Adiponectin
Reactivity:	Human, Mouse, Rat
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences