

CCL2 / MCP-1 Rabbit anti-Mouse Polyclonal (FITC) Antibody - LS-C212202 - LSBio

CatalogID:	LS-C212202
Target:	chemokine (C-C motif) ligand 2 (CCL2)
Synonyms:	CCL2 Antibody, C-C motif chemokine 2 Antibody, Chemokine (C-C motif) ligand 2 Antibody, GDCF-2 Antibody, HC11 Antibody, HSMCR30 Antibody, MCP1 Antibody, SCYA2 Antibody, Small-inducible cytokine A2 Antibody, Monocyte secretory protein JE Antibody, MCAF Antibody, MCP-1 Antibody, Monocyte chemotactic protein 1 Antibody, SMC-CF Antibody
Family / Subfamily:	Interocrine / Interocrine beta
Host	CCL2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CCL2 / MCP-1 antibody was raised against Mouse
Antigen Type:	Recombinant protein
Immunogen:	CCL2 / MCP-1 antibody was raised against recombinant mouse C-C motif chemokine 2 protein.
Specificity:	Mouse CCL2 / MCP-1
Reactivity:	Mouse
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences