

VEGFC / VEGF-C Rabbit anti-Mouse Polyclonal (FITC) Antibody - LS-C212207 - LSBio

CatalogID:	LS-C212207
Target:	vascular endothelial growth factor C (VEGFC)
Synonyms:	VEGFC Antibody, Flt4-L Antibody, FLT4 ligand DHM Antibody, VEGF-C Antibody, VRP Antibody, Flt4 ligand Antibody
Family / Subfamily:	PDGF / not assigned-PDGF
Host	VEGFC antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	VEGFC / VEGF-C antibody was raised against Mouse
Antigen Type:	Recombinant protein
Immunogen:	VEGFC / VEGF-C antibody was raised against recombinant Mouse Vascular endothelial growth factor C protein.
Specificity:	Mouse VEGFC / VEGF-C
Reactivity:	Mouse
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences