

ALB / Serum Albumin Rabbit anti-Pigeon Polyclonal (FITC) Antibody - LS-C212225 - LSBio	
CatalogID:	LS-C212225
Target:	albumin
Synonyms:	ALB Antibody, Albumin (32 AA) Antibody, Albumin (AA 34) Antibody, Albumin Antibody, Growth-inhibiting protein 20 Antibody, PRO1341 Antibody, PRO0883 Antibody, Serum albumin Antibody, PRO0903 Antibody
Host	ALB antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ALB / Serum Albumin antibody was raised against Pigeon
Antigen Type:	Purified protein
Immunogen:	ALB / Serum Albumin antibody was raised against serum Albumin purified from pigeon serum.
Specificity:	Pigeon ALB / Serum Albumin
Reactivity:	Pigeon
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	