

IL1A / IL1 Alpha Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C212324 - LSBio

CatalogID:	LS-C212324
Target:	interleukin 1, alpha (IL1A)
Synonyms:	IL1A Antibody, Hematopoietin-1 Antibody, IL-1 alpha Antibody, IL1 Antibody, IL1-ALPHA Antibody, IL1 Alpha Antibody, Interleukin-1 alpha Antibody, IL-1A Antibody, IL1F1 Antibody, Pro-interleukin-1-alpha Antibody, Preinterleukin 1 alpha Antibody, Interleukin 1, alpha Antibody
Family / Subfamily:	IL-1 / not assigned-IL-1
Host	IL1A antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	IL1A / IL1 Alpha antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	IL1A / IL1 Alpha antibody was raised against recombinant human Interleukin-1 alpha protein.
Specificity:	Human IL1A / IL1 Alpha
Reactivity:	Human
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences