

AKT2 Rabbit anti-Mouse Polyclonal (aa92-122) (FITC) Antibody - LS-C212658 - LSBio

CatalogID:	LS-C212658
Target:	v-akt murine thymoma viral oncogene homolog 2 (AKT2)
Synonyms:	AKT2 Antibody, PKB beta Antibody, PKBBETA Antibody, PRKBB Antibody, Protein kinase Akt-2 Antibody, Protein kinase B beta Antibody, RAC-PK-beta Antibody, Rac protein kinase beta Antibody, RAC-BETA Antibody, HIHGHH Antibody, PKBB Antibody
Family / Subfamily:	Protein Kinase / RAC/AKT
Host	AKT2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	AKT2 antibody was raised against Mouse
Antigen Type:	Synthetic peptide
Immunogen:	AKT2 antibody was raised against kLH-conjugated synthetic peptide (aa92-122) from the N-terminal region of mouse Akt2.
Specificity:	Mouse AKT2
Epitope:	aa92-122
Reactivity:	Mouse, Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences