

**AP1M1 Rabbit anti-Mouse Polyclonal (aa212-242) (AP) Antibody - LS-C212733 - LSBio**

<b>CatalogID:</b>	LS-C212733
<b>Target:</b>	Adaptor-related Protein Complex 1, Mu 1 Subunit (AP1M1)
<b>Synonyms:</b>	AP1M1 Antibody, AP47 Antibody, Clathrin adaptor protein AP47 Antibody, HA1 47 kDa subunit Antibody, Mu1A-adaptin Antibody, AP-1 complex subunit mu-1 Antibody, AP-mu chain family member mu1A Antibody, CLAPM2 Antibody, CLTNM Antibody, MU-1A Antibody, Mu-adaptin 1 Antibody
<b>Host</b>	AP1M1 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Alkaline Phosphatase (AP)
<b>Immunogen Species:</b>	AP1M1 antibody was raised against Mouse
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	AP1M1 antibody was raised against kLH-conjugated synthetic peptide (aa212-242) of mouse AP1M1.
<b>Specificity:</b>	Mouse AP1M1.
<b>Epitope:</b>	aa212-242
<b>Reactivity:</b>	Mouse
<b>Purification:</b>	Protein A purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences