

AP1M1 Rabbit anti-Mouse Polyclonal (aa212-242) (APC) Antibody - LS-C212734 - LSBio

CatalogID:	LS-C212734
Target:	Adaptor-related Protein Complex 1, Mu 1 Subunit (AP1M1)
Synonyms:	AP1M1 Antibody, AP47 Antibody, Clathrin adaptor protein AP47 Antibody, HA1 47 kDa subunit Antibody, Mu1A-adaptin Antibody, AP-1 complex subunit mu-1 Antibody, AP-mu chain family member mu1A Antibody, CLAPM2 Antibody, CLTNM Antibody, MU-1A Antibody, Mu-adaptin 1 Antibody
Host	AP1M1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Allophycocyanin (APC)
Immunogen Species:	AP1M1 antibody was raised against Mouse
Antigen Type:	Synthetic peptide
Immunogen:	AP1M1 antibody was raised against kLH-conjugated synthetic peptide (aa212-242) of mouse AP1M1.
Specificity:	Mouse AP1M1.
Epitope:	aa212-242
Reactivity:	Mouse
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences