

AP1M1 Rabbit anti-Human Polyclonal (aa206-235) (Azide-free) (HRP) Antibody - LS-C212743 - LSBio	
CatalogID:	LS-C212743
Target:	Adaptor-related Protein Complex 1, Mu 1 Subunit (AP1M1)
Synonyms:	AP1M1 Antibody, AP47 Antibody, Clathrin adaptor protein AP47 Antibody, HA1 47 kDa subunit Antibody, Mu1A-adaptin Antibody, AP-1 complex subunit mu-1 Antibody, AP-mu chain family member mu1A Antibody, CLAPM2 Antibody, CLTNM Antibody, MU-1A Antibody, Mu-adaptin 1 Antibody
Host	AP1M1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Horseradish Peroxidase (HRP)
Modifications:	Azide-free
Immunogen Species:	AP1M1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	AP1M1 antibody was raised against kLH-conjugated synthetic peptide (aa206-235) of human AP1M1.
Specificity:	Human AP1M1
Epitope:	aa206-235
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2. No preservative added (azide free)
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	