

MAP1LC3B / LC3B Rabbit anti-Human Polyclonal (aa8-33) (FITC) Antibody - LS-C212763 - LSBio	
CatalogID:	LS-C212763
Target:	microtubule-associated protein 1 light chain 3 beta (MAP1LC3B)
Synonyms:	MAP1LC3B Antibody, ATG8F Antibody, LC3B Antibody, MAP1A/1BLC3 Antibody, MAP1A/MAP1B LC3 B Antibody, MAP1A/MAP1B light chain 3 B Antibody, MAP1ALC3 Antibody, MAP1LC3B-a Antibody
Host	MAP1LC3B antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	MAP1LC3B / LC3B antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	MAP1LC3B / LC3B antibody was raised against kLH-conjugated synthetic peptide (aa8-33) from the N-terminal region of human APG8b(MAP1LC3B).
Specificity:	Human MAP1LC3B / LC3B
Epitope:	aa8-33
Reactivity:	Human, Mouse
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	