

CCNB1 / Cyclin B1 Rabbit anti-Human Polyclonal (pSer35) (FITC) Antibody - LS-C213120 - LSBio	
CatalogID:	LS-C213120
Target:	cyclin B1 (CCNB1)
Synonyms:	CCNB1 Antibody, CCNB Antibody, Cyclin B1 Antibody, G2/mitotic-specific cyclin-B1 Antibody, G2/mitotic-specific cyclin B1 Antibody
Host	CCNB1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CCNB1 / Cyclin B1 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	CCNB1 / Cyclin B1 antibody was raised against kLH-conjugated synthetic phosphopeptide around Ser35 of human CCNB1.
Specificity:	Human CCNB1 / Cyclin B1
Epitope:	pSer35
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	