

CD40 Rabbit anti-Human Polyclonal (aa246-276) (FITC) Antibody - LS-C213173 - LSBio

CatalogID:	LS-C213173
Target:	CD40 molecule, TNF receptor superfamily member 5
Synonyms:	CD40 Antibody, B-cell surface antigen CD40 Antibody, B cell-associated molecule Antibody, B cell surface antigen CD40 Antibody, CD40 antigen Antibody, CD40 type II isoform Antibody, p50 Antibody, Bp50 Antibody, CD40L receptor Antibody, CDW40 Antibody, TNFRSF5 Antibody
Family / Subfamily:	TNF Receptor / not assigned-TNF Receptor
Host	CD40 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD40 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CD40 antibody was raised against kLH-conjugated synthetic peptide (aa246-276) from the C-terminus of human CD40.
Specificity:	Human CD40
Epitope:	aa246-276
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences