

CD45 Rabbit anti-Human Polyclonal (aa1252-1282) (APC) Antibody - LS-C213176 - LSBio

CatalogID:	LS-C213176
Target:	protein tyrosine phosphatase, receptor type, C (PTPRC)
Synonyms:	PTPRC Antibody, B220 Antibody, CD45R Antibody, CD45 Antibody, gp180 Antibody, Leukocyte common antigen Antibody, L-CA Antibody, Leucocyte common antigen Antibody, LY5 Antibody, T200 Antibody, T200 glycoprotein Antibody, T200 leukocyte common antigen Antibody, CD45 antigen Antibody, LCA Antibody
Family / Subfamily:	Protein Phosphatase / R1/R6
Host	PTPRC antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Allophycocyanin (APC)
Immunogen Species:	CD45 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CD45 antibody was raised against kLH-conjugated synthetic peptide (aa1252-1282) from the C-terminus of human CD45.
Specificity:	Human CD45.
Epitope:	aa1252-1282
Reactivity:	Human
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.03% sodium azide.
Recommended Storage:	May be stored at 4°C before opening. DO NOT FREEZE! Stable at 4°C as an undiluted liquid. Stable for at least 12 months at 4°C. Sensitive to light.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Concentration:	0.5 mg/ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences