

**CD45 Rabbit anti-Human Polyclonal (aa1252-1282) (PE) Antibody - LS-C213199 - LSBio**

<b>CatalogID:</b>	LS-C213199
<b>Target:</b>	protein tyrosine phosphatase, receptor type, C (PTPRC)
<b>Synonyms:</b>	PTPRC Antibody, B220 Antibody, CD45R Antibody, CD45 Antibody, gp180 Antibody, Leukocyte common antigen Antibody, L-CA Antibody, Leucocyte common antigen Antibody, LY5 Antibody, T200 Antibody, T200 glycoprotein Antibody, T200 leukocyte common antigen Antibody, CD45 antigen Antibody, LCA Antibody
<b>Family / Subfamily:</b>	Protein Phosphatase / R1/R6
<b>Host</b>	PTPRC antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	CD45 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	CD45 antibody was raised against kLH-conjugated synthetic peptide (aa1252-1282) from the C-terminus of human CD45.
<b>Specificity:</b>	Human CD45.
<b>Epitope:</b>	aa1252-1282
<b>Reactivity:</b>	Human
<b>Purification:</b>	Protein G purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences