

COLQ Rabbit anti-Human Polyclonal (aa324-353) (FITC) Antibody - LS-C213380 - LSBio

CatalogID:	LS-C213380
Target:	collagen-like tail subunit (single strand of homotrimer) of asymmetric acetylcholinesterase (COLQ)
Synonyms:	COLQ Antibody, AChE Q subunit Antibody, EAD Antibody
Family / Subfamily:	Collagen / not assigned-Collagen
Host	COLQ antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	COLQ antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	COLQ antibody was raised against kLH-conjugated synthetic peptide (aa324-353) from the C-terminus of human COLQ.
Specificity:	Human COLQ
Epitope:	aa324-353
Reactivity:	Human, Mouse
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences