

COPE Rabbit anti-Human Polyclonal (aa280-308) (PE) Antibody - LS-C213388 - LSBio

CatalogID:	LS-C213388
Target:	coatomer protein complex, subunit epsilon (COPE)
Synonyms:	COPE Antibody, Coatomer subunit epsilon Antibody, Coatomer epsilon subunit Antibody, Epsilon coat protein Antibody, Epsilon-coat protein Antibody, Epsilon-COP Antibody
Host	COPE antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	COPE antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	COPE antibody was raised against kLH-conjugated synthetic peptide (aa280-308) from the C-terminus of human COPE.
Specificity:	Human COPE
Epitope:	aa280-308
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences